

Post- Use Defibrillator Checklist

***Each time the defibrillator is used, contact the Occupational Health Nurse Consultant in the OHS branch on extension 51014 immediately for replacement pads.**

Defibtech DDU-100
Defibtech DDU-100 Serial Number:
Defibtech DDU-100 Location:

Complete the following checklist:

Refer to Defibtech User Manual, sections 4.9 & 5.1-5.8 for more detailed instructions (Manual should be kept with defibrillator). Contact OHS for assistance with the below if required.

1. Unit and accessories are free from damage dirt and contamination. Clean or replace if necessary (refer to manual).	<input type="checkbox"/> Please tick when complete	Date/...../.....
2. Battery pack has not past expiration date (remove/reinsert).	<input type="checkbox"/> Please tick when complete	Date/...../.....
3. Pads replaced and within expiry date (See * above).	<input type="checkbox"/> Please tick when complete	Date/...../.....
4. Manually initiated self-test conducted.	<input type="checkbox"/> Please tick when complete	Date/...../.....
5. Active Status Indicator (ASI) is flashing green.	<input type="checkbox"/> Please tick when complete	Date/...../.....

Signature **Please print name.....**

Copy of completed form to be sent to Occupational Nurse Consultant, OHS Branch, Clayton campus –Original to be filed in area.