

MICROMON SANGER SEQUENCING

MODIFIED ALKALINE LYSIS/PEG TEMPLATE PREPARATION PROTOCOL

Reagents required:

Chloroform

Deionized water

Ethanol, 70%

GET buffer (50 mM glucose, 10 mM EDTA, 25 mM Tris, pH 8.0)

Isopropanol, 100% (anhydrous)

PEG 8000, 13% (sterilized by autoclaving, rather than by filtration)

Potassium acetate, 3 M, pH 4.8

RNase A (DNase-free), 10 mg/ml

Sodium chloride (NaCl), 4 M

Sodium hydroxide (NaOH), 0.2 N, with 1% SDS (freshly made)

To minimize shearing of contaminating chromosomal DNA, do not use a vortex during this procedure.

1 Pellet 1.5-ml aliquots of culture for one minute in a microcentrifuge at maximum speed.

Note A total culture volume of 4.5 ml can be spun down per tube without changing volumes in the procedure. This allows you to achieve a threefold increase in yield while eliminating the need for extra tubes and additional handling.

2 Remove the supernatant by aspiration.

3 Resuspend the bacterial pellet in 200 μ l of GET buffer by pipetting up and down.

4 Add 300 μ l of freshly prepared 0.2 N NaOH/1% SDS. Mix the contents of the tube by inversion. Incubate on ice for 5 minutes.

5 Neutralize the solution by adding 300 μ l of 3.0 M potassium acetate, pH 4.8. Mix by inverting the tube. Incubate on ice for 5 minutes.

6 Remove cellular debris by spinning in a microcentrifuge at maximum speed for 10 minutes at room temperature. Transfer the supernatant to a clean tube.

7 Add RNase A (DNase-free) to a final concentration of 20 μ g/ml. Incubate the tube at 37 °C for 20 minutes.

8 Extract the supernatant twice with chloroform:

- a. Add 400 μ l of chloroform.
- b. Mix the layers by inversion for 30 seconds.
- c. Centrifuge the tube for 1 minute to separate the phases.
- d. Transfer the upper aqueous phase to a clean tube.

9 Add an equal volume of 100% isopropanol. Mix the contents of the tube by inversion.

- 10** Spin the tube in a microcentrifuge at maximum speed for 10 minutes at room temperature.
- 11** Remove the isopropanol completely by aspiration.
- 12** Wash the DNA pellet with 500 μ l of 70% ethanol. Dry under vacuum for 2–3 minutes.
- 13** Dissolve the pellet in 32 μ l of deionized water.
- 14** Add 8.0 μ l of 4 M NaCl, then 40 μ l of autoclaved 13% PEG 8000.
- 15** Mix thoroughly, then leave the sample on ice for 20 minutes.
- 16** Pellet the plasmid DNA by spinning in a microcentrifuge at maximum speed for 15 minutes at 4–8 °C.
- 17** Carefully remove the supernatant. Rinse the pellet with 500 μ l of 70% ethanol.
- 18** Resuspend the pellet in 20 μ l of deionized water. Store at -20 °C.